

BOECO VACUUM PUMP R-300

is a piston-powered, oil-free vacuum pump, adaptable to many laboratory requirements.

Quiet and low vibration

Directly motor driven without transmission mechanism and quality rubber feet keep the noise level only at approx 54 dB.

Compact and light weight.

Aluminium die-cast and precise construction make the R-300 so small and compact, weight only 4,1 kg

Clean and maintenance free

Oil-free design makes the R-300 clean and maintenance free, we guarantee two years of 3000 working hours of free service parts (Excluding moisture filter)



Code	Description
BOE 8830000	BOECO R-300 vacuum pump 220V, 50/60 Hz

Specification	
Max. vacuum:	600 mm hg
Vacuum speed (2,5 L tank):	50 Hz
200 mmHg:	3 sec
300 mmHg:	5 sec
400 mmHg:	8 sec
500 mmHg:	14 sec
600 mmHg:	25 sec
Max flow rate:	13 l/min
Motor rotation:	1400 rpm
Horse power:	1/8 HP
Pole:	4 P
Net weight:	4.1.kgs
Port thread:	5/16 inch
Noise level:	54db